

NexTemp® Ultra

CLINICALLY ACCURATE, SINGLE-USE THERMOMETERS OFFERING
UNSURPASSED ACCURACY & MAXIMUM INFECTION CONTROL

PRODUCT SPECIFICATIONS

- Clinical accuracy to ± 0.1 °C/ ± 0.2 °F
- Offers a 5-year shelf life
- No calibrating or cleaning required
- Light-weight and compact
- No temperature-controlled storage needed
- Available in both Celsius and Fahrenheit
- Does not require batteries, cords, probe covers, or cleaning supplies
- Non-toxic; 100% hygienic; hypoallergenic; glass and mercury- free; not made with natural rubber latex; waterproof; unbreakable
- Featuring an extended temperature read time of 60 seconds

PACKAGING OPTIONS

NexTemp® Ultra single-use thermometers are available in Clinical Boxes of 100, as well as Individual Packs of 12 or 24 thermometers.

Our boxes of 100 are ideal for use in hospitals, clinics, urgent care centers, and teaching institutions. They're also great for large-scale temperature monitoring needs, such as return-to-work programs, conferences and expos, or sporting and entertainment events.

Our Individual Packages come in either 12 or 24 thermometers per pack. And we offer both stock designs, as well as the option to fully customize the packaging! You can add your logo, colors, wording, and imagery, and they can be produced in any language for international use! These packs are ideal for at-home use, in first-aid or medical supply kits, travel programs, disaster relief initiatives, and are widely sold in retail settings, both in-store and on-line.

THE BENEFITS OF DISPOSABLES

When thermometers are used on multiple patients, the potential for transmission of germs and pathogens significantly increases. NexTemp® Ultra single-use thermometers are individually wrapped and intended to be disposed of after use, circumventing the risk of infection transfer by providing a clean instrument that is used once and then discarded.

