

Tempa•DOT®

**CLINICALLY ACCURATE, SINGLE-USE THERMOMETERS OFFERING
 UNSURPASSED ACCURACY & MAXIMUM INFECTION CONTROL**

PRODUCT SPECIFICATIONS

- Clinical accuracy to ± 0.1 °C/ ± 0.2 °F
- Offers a 2-year shelf life
- No calibrating or cleaning required
- Light-weight and compact
- Should be stored at or below 30°C/86°F
- Available in both Celsius and Fahrenheit
- Does not require batteries, cords, probe covers, or cleaning supplies
- Non-toxic; 100% hygienic; hypoallergenic; glass and mercury- free; not made with natural rubber latex; waterproof; unbreakable
- Temperature reading is permanently locked in
- Part No: 5532

ABOUT TEMPADOT® THERMOMETERS

The Tempa•DOT® line of clinically accurate, single-use thermometers was developed to answer the need for a thermometer that provides highly accurate readings, while also protecting patients and healthcare providers from the spread of germs and infections. All of our Tempa•DOT® thermometers are 100% hygienic and hypoallergenic. And they are manufactured in our FDA and ISO Certified facilities, in a controlled and regulated environment, to ensure their safety.

ADDED PROTECTION WITH STERILIZATION

In addition to our standard Tempa•DOT® thermometers, we also offer Tempa•DOT® Sterile. Tempa•DOT® Sterile thermometers are identical in design and functionality as our standard thermometers that have also been sterilized for instances where infection control is of critical importance. Tempa•DOT® Sterile is widely used in hospitals and medical facilities around the world, as well as dialysis centers, chemotherapy infusion clinics, and ICU wards, among others.

