

# Tempa•DOT®

CLINICALLY ACCURATE, SINGLE-USE THERMOMETERS OFFERING  
UNSURPASSED ACCURACY & MAXIMUM INFECTION CONTROL

## PRODUCT SPECIFICATIONS

- Clinical accuracy to  $\pm 0.1$  °C/ $\pm 0.2$  °F
- Offers a 2-year shelf life
- No calibrating or cleaning required
- Light-weight and compact
- Should be stored at or below 30°C/86°F
- Available in both Celsius and Fahrenheit
- Does not require batteries, cords, probe covers, or cleaning supplies
- Non-toxic; 100% hygienic; hypoallergenic; glass and mercury- free; not made with natural rubber latex; waterproof; unbreakable
- Temperature reading is permanently locked in
- Part No: 5532NS

## ABOUT TEMPADOT® THERMOMETERS

The Tempa•DOT® line of clinically accurate, single-use thermometers was developed to answer the need for a thermometer that provides highly accurate readings, while also protecting patients and healthcare providers from the spread of germs and infections. All of our Tempa•DOT® thermometers are 100% hygienic and hypoallergenic. And they are manufactured in our FDA and ISO Certified facilities, in a controlled and regulated environment, to ensure their safety.

## ENSURING SAFETY WITH BIOBURDEN TESTING

Safety is paramount at Medical Indicators. To ensure their safety, we test our 5532NS Tempa•DOT® thermometers utilizing bioburden testing. The bioburden rate on our 5532NS thermometers is <1CFU - meaning that our thermometers have a population of less than one viable microorganism and show no detectable bacterial or fungal bioburden.

